

Cambridgeshire and Peterborough Wheelchair Service

Covid-19 – Cleaning of Wheelchairs

This document is intended for use by service users, carers, family members, PAs or anyone responsible for a wheelchair user/the maintenance of their wheelchair.

Regular cleaning of your wheelchair should already be in place, but there are a few extra steps you can take during the Covid-19 pandemic.

These are:

- Use wipes with an alcohol content of at least 70%
- Clean the wheels, arm rests and controller before entering the house
- Wash any removable covers and cushions as per manufacturer instructions
- If you can, store the wheelchair outside of the house (including garage, shed or other outhouse if it is raining)
- Before using the chair again, wipe down all touchable surfaces once again

Service provided by

